

Human Mesenchymal Stem Cell Adipogenesis Detection Kit (HMSC-A PCR)

Cat. No. 8298, 50 reactions

Product Description

Human mesenchymal stem cells (MSC) are a population of multipotent cells that can be differentiated into multiple lineage-specific cells, which can form bone, fat, cartilage, muscle and tendon. Among them, adipocytes are one of the cell types that can be derived from MSCs in the process of adipogenesis. The process of fat formation plays a role in obesity, cardiovascular diseases and metabolic disorders.

ScienCell has created a convenient multiplex PCR kit for the routine detection of human mesenchymal stem cell adipogenesis. Multiplex PCR allows two or more primers to be amplified in a single PCR reaction by using multiple primer pairs in a single reaction mixture, allowing for considerable saving in labor, cost and precious DNA samples. All required PCR reagents are supplied in this kit. Simply add DNA template and perform PCR reaction. Tube 1 ready-mix reaction contains *PPARy* and *FABP4* primers that allow for the detection of mid- and late-stage of human MSC adipogenesis, respective.

Kit Components

| Cat. No. | # of vials | Name | Quantity | Storage |
|----------|------------|-----------------------------------|----------|---------|
| 8298a | 1 | Tube 1: differentiation ready-mix | 660 μL | -20°C |
| 8298b | 1 | Nuclease-free H ₂ O | 1 mL | -20°C |

Materials to be Supplied by the User

DNA templates Thin wall PCR tubes Thermal cycler Agarose gel Ethidium bromide Electrophoresis system Gel imager

Quality Control

cDNAs from differentiated human mesenchymal stem cells during the time-course of adipogenesis were used as template DNA. Each PCR product was sequenced to ensure specificity.

Product Use

HMSC-A PCR kit is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store at -20°C upon receipt. Avoid repeated freeze thaw cycles by making six aliquots at 110 µL each.

Shipping

Dry ice.

1. Mix the following components in a thin-wall PCR tube:

| Component | Amount | |
|--------------------------------|-------------|--|
| Tube 1 ready-mix | 12 μL | |
| cDNA template* | 2 μg | |
| Nuclease-free H ₂ O | up to 20 μL | |

^{*}cDNA concentration should be ~2 μg/μL

2. Perform PCR using the following conditions:

 Step 1: 95°C
 3 min

 Step 2: 95°C
 30 sec


 Step 3: 58°C
 30 sec

 Step 4: 72°C
 1 min

Step 5: Repeat Steps 2-4 for 24 times

Step 6: 72°C 2 min

3. Visualize PCR products on a 1.5% agarose gel containing ethidium bromide.



4. Expected product sizes:

| Gene | Expected Size | | |
|--------------|---------------|--|--|
| $PPAR\gamma$ | 505bp | | |
| FABP4 | 251bp | | |