

# Human Pulmonary Artery Endothelial Cell Total RNA (HPAEC tRNA)

Catalog #3105

Quantity: 10 µg

# **Product Description**

Human Pulmonary Artery Endothelial Cell total RNA (HPAEC tRNA) is prepared from early passage Human Pulmonary Artery Endothelial Cells using the Qiagen AllPrep DNA/RNA Mini kit. Total RNA concentration was determined by DU 800 UV/Visible spectrophotometer. Total RNA from ScienCell Research Laboratories is convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues for isolation of RNA.

# **Quality Control**

- Purity: A<sub>260</sub>/A<sub>280</sub> ratio of 1.9-2.1
- RNA integrity: intact 18S rRNA bands and 28S rRNA bands verified by a denaturing agarose gel, with a 2:1 ratio of 28S:18S.

## **Application**

- RT-PCR and Real Time RT-PCR
- cDNA synthesis
- RNA protection assay
- Northern blot analysis
- Purification of mRNA for library construction
- PCR-select cDNA subtraction

#### **Product Use**

HPAEC tRNA is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

#### **Storage**

Store the RNA immediately at -80°C upon arrival. To ensure stability, it is recommended to aliquot the RNA into single-use fractions to avoid repeated freeze/thaw cycles.

### Shipping

Dry ice.

Note: The Qiagen AllPrep DNA/RNA Mini kit isolates all RNA molecules larger than 200 nucleotides. RNA molecules <200 nucleotides are selectively excluded.