

# Human Ovarian Surface Epithelial Cell Lysate (HOSEpiC Lysate) Catalog #7316

Quantity: 200 µg

### **Product Description**

Human Ovarian Surface Epithelial Cell Lysate (HOSEpiC Lysate) is prepared from early passage Human Ovarian Surface Epithelial Cells using a modified RIPA buffer (50 mM Tris-HCl pH 7.2, 150 mM NaCl, 1% NP-40, 1mM EDTA, 1 mM EGTA, 0.4 mM PMSF, 5  $\mu$ g/ml aprotinin, 5  $\mu$ g/ml leupeptin, 1  $\mu$ g/ml pepstatin, 1 mM Na<sub>3</sub>VO<sub>4</sub>, 5 mM NaF). Lysate concentration is determined by Beckman Coulter DU 800 spectrophotometer. Lysate is diluted and stored in Storage Buffer (Tris-HCl pH 6.8, 5% glycerol, 1% SDS, 0.002% bromphenol blue, 5%  $\beta$ -mercaptoethanol). Lysate from ScienCell Research Laboratories is convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

### **Quality Control**

• The quality of the lysate is confirmed by SDS-PAGE.

### Application

• Western Blotting

## **Product Use**

HOSEpiC Lysate is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### Storage

Store lysate immediately at -80°C upon arrival; aliquot the lysate to avoid repeated freeze/thaw cycles

### Shipping

Dry ice.

Note: It is recommended to heat the lysate for 5 min at 95°C and put on ice before loading to a SDS-PAGE gel.